


Přehled technických novinek z pohledu společnosti Siemens



The Biograph family

Innovation leadership in PET/CT



Some of these products and features are not commercially available in all countries. Future availability cannot be guaranteed.

¹ Fastest measured value on a single system. Based on competitive literature available at time of publication. Data on file.

The Biograph family

Our PET/CT scanner portfolio


PMT

Biograph™ Horizon



More within reach

Biograph Vision™



See a whole new world of precision.



Biograph mCT



A wide bore and more

Biograph Vision Quadra™



Bigger perspective. Better answers.

These products and features are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Biograph Vision PET/CT family



Biograph Vision™ 450



Biograph Vision 600



Biograph Vision.X™



Biograph Vision and Biograph Vision.X are not commercially available in all countries. Future availability cannot be guaranteed.

Unrestricted © Siemens Healthcare GmbH, 2018 MI-4180

Biograph Vision

Going beyond digital with the Optiso UDR detector



3.2 mm
LSO
crystals

214/178
ps¹
time of
flight

100 cps/kBq¹
effective
sensitivity


100%
SiPM
sensor
coverage¹

Accuracy
Performance
Reproducibility

¹ Based on internal measurements at time of publication. Data on file;
Biograph Vision and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed.
Please contact your local Siemens Healthineers organization for further details.

Accuracy

Accurately identifying the appropriate disease stage is key to defining the right treatment



Disease location

Local?

Regional?

Distant?

Therapy option

Surgery?

Radiation?


Chemo?

Today's challenge

- Identifying additional local lesions, lymph node involvement or small metastases can change patient staging and influence therapy decisions
- **Poor spatial resolution can affect**
 - Lesion detectability
 - Image quality

Insufficient spatial resolution can negatively affect lesion detectability and staging.

Transcend digital with Optiso UDR detector



LSO
scintillator



3.2-mm
crystals



SiPM technology with
100% coverage¹



Optiso UDR
detector

¹ Based on internal measurements at time of publication. Data on file.

Optiso

UDR detector

revolutionizing precision in PET/CT

We are determined to do more than just continue the momentum of improving technology in increments. We're inspired to make breakthrough improvements that would revolutionize PET/CT. The result?

Transcending digital with the Optiso Ultra Dynamic Range (UDR) detector.

- Optiso UDR detector technology reveals a new world of precision—to help you detect small lesions and devise accurate treatment strategies


1

2

3

4

5



1. Lutetium oxyorthosilicate (LSO), a fast and efficient scintillator, is grown and cut in-house through a vertically integrated manufacturing process to ensure the highest quality.
2. 3.2-mm crystal elements are individually selected and deliver high isotropic spatial resolution; higher spatial resolution may result in improved lesion detectability.
3. 100% coverage¹ of the crystal area with SiPM sensors results in a timing resolution of 214-ps¹ and 3.9x higher effective sensitivity² for faster scans and lower dose
4. A small block size delivers higher¹ effective peak NECR for improved clinical performance.²
5. High-flow direct-cooling of the detector plate allows the detector to operate at room temperature¹ for outstanding performance, serviceability and improved patient comfort.

¹ Based on internal measurements at time of publication. Data on file.

² Compared to current Siemens state-of-the art technologies. Data on file.


Reducing crystal size improves spatial resolution and detectability

3.2-mm crystals


for improved, more precise imaging

Smaller crystal elements and block size improve small-lesion detectability by delivering better volumetric resolution. 3.2-mm crystals allow you to see smaller lesions to help you confidently stage, risk-stratify, and develop appropriate treatment strategies sooner.

- Small crystal elements improve signal-to-noise ratio and quantitative accuracy
- 60% better volumetric resolution¹




60%
Better volumetric resolution¹



¹ Compared to competitive literature available at time of publication. Data on file.

Clear visualization of the 2.4 mm cylinders of the Mini-Derenzo phantom with Biograph Vision




SiPM technology with 100% coverage

for optimized performance

Time-of-flight (TOF) performance depends on collecting light from all photons in the scintillation. The Optiso UDR is designed so SiPMs cover the entire lutetium oxyorthosilicate (LSO)-array area, allowing all light from the scintillation to be detected. This leads to 100% coverage¹ and enables fast temporal resolution.

- Biograph Vision™ provides 214-ps (picosecond) temporal resolution¹ for best-in-class TOF and effective sensitivity²
- Biograph Vision 3.9x TOF gain³ amplifies scanner sensitivity for faster scans and lower dose



¹ Based on internal measurements at time of publication. Data on file.

² Compared to competitive literature available at time of publication. Data on file

³ Compared to current Siemens Healthineers state-of-the-art technologies. Data on file.



Biograph Vision is not commercially available in all countries. Its future availability cannot be guaranteed.

Time of flight (TOF) PET

TOF is the ability to measure the time difference between the arrival of the two coincidence photons in the PET detector and use that information to:

- Improve lesion detectability and overall image quality
- Improve clinical efficiency by enhancing sensitivity, enabling faster scans and lower dose
- Improve effectiveness of all PET corrections (attenuation, scatter, randoms, MoCo, etc) and image reconstruction

Faster time of flight delivers higher sensitivity gain




Direct histogram PET image
Image reconstruction input



TOF gain
D = 20 cm phantom diameter

$$TOF_{Gain} = \frac{2D}{c\Delta t}$$


D = object size
c = speed of light
Delta (t) = TOF resolution




Budinger, Thomas F., "Time-of-Flight Positron Emission Tomography: Status Relative to Conventional PET", *Journal of Nuclear Medicine*, 1983;24:73-78. 14

Faster TOF reduces noise, which increases detectability of small lesions


non TOF



TOF (540ps)



TOF (214ps)



High-resolution torso phantom

Sphere size (mm): 5.0, 7.9, 9.9, 12.4, 15.4, 19.8

6:1 contrast-to-background

Left and middle image acquired on Biograph mCT. Right image acquired on Biograph Vision.

Biograph Vision and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed.

Please contact your local Siemens Healthineers organization for further details.

Performance

Image quality or patient dose?

An ongoing compromise



Today's challenge

- Delivering high diagnostic value images at lower doses
- Addressing throughput demands

Balancing speed, dose and image quality requires the user to compromise somewhere.

Biograph Vision's high performance enables high-quality scans with less dose and in less time

3.9X

higher effective
sensitivity¹

"We have now already, as compared to the older system, reduced the activity we inject... Now it's probably 30% faster with about 30% less dose which is something very acceptable."¹

Prof. John Prior, MD, PhD

Head of Department, Nuclear Medicine

Centre Hospitalier Universitaire Vaudois CHUV, Lausanne, Switzerland

Significantly reduce:

- Injected dose
- Scan time

without compromising image quality

¹ Compared to current state-of-the-art technologies. Data on file.



The statements by Siemens' Healthineers customers described herein are based on results that were achieved in the customer's unique setting. Because there is no "typical" hospital or laboratory and many variables exist (e.g., hospital size, samples mix, case mix, level of IT and/or automation adoption) there can be no guarantee that other customers will achieve the same results.

Biograph Vision and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed.

Please contact your local Siemens Healthineers organization for further details.

Faster time-of-flight amplifies scanner sensitivity for faster scans and lower doses

Effective sensitivity = NEMA Sensitivity * TOF Gain



¹Based on internal measurements available at time of publication. Data on file. ²Gain calculated for a 20 cm cylindrical object. Biograph Vision and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.


Better effective sensitivity improves total cost of ownership by enabling scans with reduced tracer dose

43% reduction of FDG injected dose compared to previous system at Centre Hospitalier Universitaire Vaudois (CHUV)

Changed the protocol from 3.5 MBq/kg (0.043 mCi/lb) to 2 MBq/kg (0.025 mCi/lb)

4.4 mCi for a 180 lb patient

Dose cost savings*
(accumulated in thousand \$)



* Savings simulated assuming a cost of \$ 200.00 per dose and an average of 50 patients/week

The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

Biograph Vision and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Reproducibility



Low resolution can prevent new lesions from being detected and can lead to inaccurate quantitative values

146k

Yearly cost of treatment with Imatinib mesylate (IM) in US dollars¹

15%

of patients displayed primary resistance to IM treatment²

Discerning non-responders as early as possible in cancer treatment can maximize effective patient care and outcomes

Today's challenge

- Low resolution can impact lesion detectability
- Quantitative values are not always accurate or reproducible
- Low resolution and quantitative accuracy limit the ability to evaluate the progression of the disease and adjust treatment

¹ <http://www.ascopost.com/issues/may-25-2016/the-arrival-of-generic-imatinib-into-the-us-market-an-educational-event/> accessed on 02/Apr/2018; ²Farag, et al, Early Evaluation of Response Using 18F-FDG PET Influences Management in Gastrointestinal Stromal Tumor Patients Treated with Neoadjuvant Imatinib, J Nucl Med 2018
Biograph Vision and its features and applications are not commercially available in all countries. Their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

AIDAN Platform

Intelligent imaging platform for PET/CT

AIDAN is powered by unique technologies

ALPHA (Anatomical Landmarking and Parsing of Human Anatomy) a sophisticated AI engine supporting optimized workflows on both scanners and reading solutions.




SMART patient handling system (PHS) an exclusive bed design with magnetically drive precision. Zero differential deflection for accurate attenuation correction and TG-66 compliant to support positioning of radiation oncology devices.

AI-based features at the scanner are based on unique ALPHA technology

ALPHA uses landmarks and reference regions to learn anatomical locations and can recognize features like human visual recognition.

- Algorithm learns from a library of thousands of expert annotated examples
- Use of redundancy: large number of points, spatial relationships, and multiple scales
- Robust regardless of disease, body position and habitus, image range and quality


ALPHA is the foundation for AI-based PET applications at the scanner.



The Biograph family—with AIDAN PET/CT AI-powered applications


The AIDAN platform unlocks a host of advanced features such as:

FlowMotion™ AI¹



Standardize protocols and personalize scans through a dedicated AI algorithm that automatically defines ranges based on each patient's unique anatomy.


OncoFreeze™ AI¹



Deliver images free of motion without extended scan time or external devices—with the click of a button.


¹ These products and features are not all available on Biograph Vision Quadra. Their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details. AIDAN is not commercially available in all countries. Its future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

FlowMotion AI: Moving the standard to personalized care



The features herein are pending 510(k) clearance, and are not yet commercially available in the United States or other countries worldwide on Biograph Horizon, Biograph mCT and Biograph Vision. Their future availability cannot be guaranteed. The features herein do not yet fulfill all the essential requirements according to the European Medical Device Directive (93/42/EEC) and its national implementations for Biograph mCT and Biograph Vision. They are not yet commercially available in the EU. Future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details. Unrestricted © Siemens Healthcare GmbH, 2018 MI-4180

FlowMotion AI automatically defines scan ranges based on patient anatomy, with the click of a button



The dedicated AI algorithm, based on ALPHA technology, is designed to automatically define scan ranges based on patient anatomy, independent of the user.

- **Faster workflow¹**
- **More personalized scan¹**
- **Reproducible results from any operator**
- **Optimum range for scan settings help avoid cut-offs or over radiation**



¹ Compared to PET/CT systems without FlowMotion™

Dedicated AI algorithm enabling FlowMotion AI Based on ALPHA Technology



What ALPHA does

- **A**utomatic **L**andmarking and **P**arsing of **H**uman **A**natomy
- Derives anatomical landmarks based on CT topogram images
- Robust regardless of disease, body position and habitus, image range and quality

How it works

- Algorithm learns from a library of thousands of expert annotated examples
- Use of redundancy: Large number of points, spatial relationships, multiple scales

What it does in FlowMotion AI

Recognizes the patient anatomy and automatically defines scan ranges based on individual organ, aiming for a faster and more personalized workflow at the scanner¹




¹ Compared to PET/CT systems without FlowMotion

Data courtesy of Praxis für Fusionierte Bildgebung, Halle (Saale), Germany.

The features herein are pending 510(k) clearance, and are not yet commercially available in the United States or other countries worldwide on Biograph Horizon, Biograph mCT and Biograph Vision. Their future availability cannot be guaranteed. The features herein do not yet fulfill all the essential requirements according to the European Medical Device Directive (93/42/EEC) and its national implementations for Biograph mCT and Biograph Vision. They are not yet commercially available in the EU. Future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

OncoFreeze AI: Deviceless solution for images virtually free of motion without extending scan time




OncoFreeze AI offers the ability to provide motion-free PET images by enabling a simple check box.




The features herein are pending 510(k) clearance, and are not yet commercially available in the United States or other countries worldwide on Biograph Horizon, Biograph mCT and Biograph Vision. Their future availability cannot be guaranteed. The features herein do not yet fulfill all the essential requirements according to the European Medical Device Directive (93/42/EEC) and its national implementations for Biograph mCT and Biograph Vision. They are not yet commercially available in the EU. Future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Today's challenges to provide respiratory motion management for every patient



Number of patients affected

90% of oncological disease is located in areas subject to respiratory motion¹



Clinical Impact

Without respiratory gating, 40% of lung lesions may even go undetected²



Setup with conventional methods

Up to 11 minutes added to every respiratory gated scan



Additional scan time required

Up to 3 times longer to scan areas that can be impacted by respiratory motion



¹BIO-TECH SYSTEMS, INC. Report 2008.

² Garcia Vicente AM, et al. (18) F-FDG PET-CT respiratory gating in characterization of pulmonary lesions: approximation towards clinical indications. Ann Nucl Med. 2010 April 24 (3) 207-14

Deviceless waveform

Artificial intelligence enables the system to identify a waveform generated from respiratory motion without the need for an external device.

This eliminates the need for additional setup to acquire a respiratory waveform to generate motion-frozen or moving PET reconstructed images.



The features herein are pending 510(k) clearance, and are not yet commercially available in the United States or other countries worldwide on Biograph Horizon, Biograph mCT and Biograph Vision. Their future availability cannot be guaranteed. The features herein do not yet fulfill all the essential requirements according to the European Medical Device Directive (93/42/EEC) and its national implementations for Biograph mCT and Biograph Vision. They are not yet commercially available in the EU. Future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Dedicated AI algorithm enabling OncoFreeze AI Based on ALPHA Technology



What ALPHA does

- Automatic Landmarking and Parsing of Human Anatomy
- Derives anatomical landmarks based on CT topogram images
- Robust regardless of disease, body position and habitus, image range and quality

How it works

- Algorithm learns from a library of thousands of expert annotated examples
- Use of redundancy: Large number of points, spatial relationships, multiple scales

What does it do in OncoFreeze AI

ALPHA technology will identify specific anatomical landmarks, to determine the zone to be corrected for respiratory motion automatically for each patient.



Data courtesy of Praxis für Fusionierte Bildgebung, Halle (Saale), Germany.

The features herein are pending 510(k) clearance, and are not yet commercially available in the United States or other countries worldwide on Biograph Horizon, Biograph mCT and Biograph Vision. Their future availability cannot be guaranteed. The features herein do not yet fulfill all the essential requirements according to the European Medical Device Directive (93/42/EEC) and its national implementations for Biograph mCT and Biograph Vision. They are not yet commercially available in the EU. Future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

NEW

Introducing Biograph Vision.X

with next-level performance

SIEMENS
Healthineers



Biograph Vision.X™ is not commercially available in all countries. Future availability cannot be guaranteed.

Unrestricted © Siemens Healthineers, 2023 MI-6100

Introducing Biograph Vision.X PET/CT

Next-level PET/CT performance

NEW



TOF
178 ps¹

Biograph Vision.X™ with exclusive, high-performance detector technology

- Industry-leading **178-ps¹** time of flight (TOF)
- Outstanding **performance gain up to 20%²**
- The fastest time of flight (TOF) performance in the industry¹

¹ Fastest measured value on a single system. Based on competitive literature available at time of publication. Data on file.


² Measured value on a single system. Data on file.

Biograph Vision.X is not commercially available in all countries. Future availability cannot be guaranteed.

Faster time of flight (TOF) amplifies scanner sensitivity for faster scans and lower doses



178-ps
timing resolution¹



¹ Fastest measured value on a single system. Based on competitive literature available at time of publication. Data on file. Biograph Vision.X™ is not commercially available in all countries. Future availability cannot be guaranteed.



The Symbia family

Setting the standard for SPECT and SPECT/CT imaging



Innovative SPECT and SPECT/CT solutions

Better outcomes and lower costs




1960s
First commercial Anger camera
To image organ functionality




1970s
First dual detector SPECT
To scan twice as fast



1989
First AUTOFORM™ collimator technology
Up to 26% higher sensitivity



1995
First automatic SPECT body contour
To get the best image resolution



2003
First 3D iterative reconstruction
Fast acquisition, lower dose



2004
First integrated ECG
To optimize workflow



First diagnostic SPECT/CT
Greater anatomical detail
Symbia™ SPECT/CT



2006
First and only automated collimator changer
To optimize workflow




First and only automated quality control
To optimize workflow




2008
First and only ultra-fast cardiac SPECT
4x more counts to double the speed



2013
First and only xSPECT™ system
Integrated SPECT and CT data
First and only high-resolution bone imaging with
Accurate and reproducible quantification



2014
Smallest SPECT room size in its class
Up to 29% smaller




2016
Standardized quantification
For multiple isotopes

2017
First SPECT/CT metal artifact reduction
Image every body¹
Symbia Intevo Bold™



2020
Standardized quantification
Low-, medium-, and high-energies



2022
Symbia Pro.specta™ with myExam Companion™
Fast, 64-slice CT imaging with automatic SPECT motion correction



¹ For patients up to 227 kg (500 lb). | All claims based on competitive literature available at time of publication. Data on file. Symbia Pro.specta and xSPECT Quant ¹³¹I are not commercially available in all countries. Future availability cannot be guaranteed.

The Symbia family

Flexibility to meet your needs today and in the future

Symbia Evo™ Excel



Symbia Evo



Symbia Pro.specta™ with myExam Companion™



Core SPECT technology: HD detectors, AUTOFORM™ collimators, high sensitivity and fine reconstructed resolution, pallet deflection reduction, 227 kg (500 lb) support, Autocontour, detector flexibility

Clinical productivity: Automatic Collimator Changer, Automatic Quality Control, IQ•SPECT™

Automatic SPECT motion
correction


Low-dose diagnostic
CT: 32/64 with iMAR

Intelligent imaging with
myExam Companion

xSPECT™ technologies:
xSPECT Bone™, xSPECT Quant™, Broad Quantification™

Symbia Pro.specta SPECT/CT with myExam Companion

Intelligent imaging in SPECT/CT



Symbia Pro.specta™ and its features are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Spending hours per week handling and moving equipment negatively impacts workplace satisfaction

Top three priorities

of nuclear medicine departments

- **Improve patient satisfaction**
- **Improve workflow and productivity**
- **Manage department costs**

Time unnecessarily spent with equipment takes technologists away from patients.

Symbia Productivity Package

Optimize workflow. Increase patient satisfaction.

Symbia™ Productivity Package

streamlines your workflow by automating collimator change and quality control, so more time can be spent performing patient-centric activities.



Automatic Collimator Changer

- Operator independent
- Less risk of equipment damage
- Conserve space
- 5 minutes gained per exchange



Automatic Quality Control

- Quality control runs overnight
- Results in the morning
- Reduced operator exposure
- Patient-focused start to each day

Automatic Quality Control

Gain up to one hour each day for more value-added tasks

Convert hundreds of hours each year from manually performing quality control to helping your patients.

Automatic Quality Control (AQC)

- Automatically run quality control overnight
- Get results first thing in the morning
- Start each day focused on your patients
- Reduce operator radiation exposure



Conventional calibration and system quality control consumes up to 1 hour per day

Daily quality control

Prepare point source	15 min
Remove collimators	05 min
Peak	01 min
Low count	12 min
Insert collimators	05 min

Total daily time spent = 38 min

Weekly quality control

Prepare point source	05 min
Remove collimators	05 min
Tuning	20 min
Insert collimators	05 min
Multi-head registration	05 min

Total weekly time spent = 40 min

Monthly quality control

Prepare point source	30 min
Multi-head registration	30 min
Remove collimators	05 min
Tuning	01 min
High count flood	120 min
Insert collimators	05 min

Total monthly time spent = 191 min

Daily + weekly + monthly

Annual workdays



1 hour
per day

Automatic Collimator Changer Simplify workflow. Eliminate variability.

Change collimators and reconfigure detectors with a just a few clicks to streamline workflow and improve the experience.

Automatic Collimator Changer (ACC)


- Change collimators up to 50% faster with ACC¹
- Eliminate variability from manual handling
- Lower risk of equipment damage
- Conserve space by reducing the need for additional collimator carts




¹Based on internal testing. Data on file.

Minimum dose. Maximum speed. The fundamentals of IQ•SPECT technology


Unique magnifying collimators



Cardio-centric acquisition






Advanced reconstruction




Advanced IQ•SPECT reconstruction

Dose protocol and acquisition time by patient

	Conventional		IQ•SPECT™	
	Dose	Time	Dose	Time
 Standard	Full	16 min.	Full	4 min.
 Pediatric	Half	16 min.	Half	8 min.
 Bariatric	Full	16 min.	Full	8 min.

Enhance your nuclear medicine department with intelligent oncology SPECT/CT imaging




¹ Compared to conventional iterative reconstruction methods

² iMAR is designed to yield images with a reduced level of metal artifacts compared to conventional reconstruction if the underlying CT data is distorted by metal being present in the scanned object. The exact amount of metal artifact reduction and the corresponding improvement in image quality achievable depends on a number of factors, including composition and size of the metal part within the object, the patient size, anatomical location and clinical practice. It is recommended, to perform iMAR reconstruction in addition to conventional reconstruction. Symbia Pro.specta™ and its features are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Symbia Pro.specta SPECT/CT with myExam Companion

MC Pro: Automatic motion correction for SPECT applications




Eliminate user variability when performing motion correction with MC Pro

- Apply motion correction automatically for SPECT applications with a simple on/off selection
- Improve image quality for most studies by automatically correcting for patient movement
- Achieve greater consistency across users and exams with automatic motion correction

Applicable for all SPECT reconstructions except for dynamic SPECT.

Symbia Pro.specta™ and its features are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

xSPECT Bone differentiates tissue boundaries for high-resolution bone imaging




xSPECT Bone™


- Uses CT as the frame-of-reference for image reconstruction
- Extracts a zone map with definition around different tissue segments to better delineate the boundaries of pharmaceutical uptake

The result: sharper definition of bony margins to generate the degree of image quality needed to reveal a patient's true condition in a single exam


xSPECT Bone image alone shows clear vertebrae uptake



CT




3D iterative with CTAC



xSPECT Bone™/CT

Reduce metal artifacts in clinical routine iMAR¹ (iterative Metal Artifact Reduction)




Reduce metal artifacts in clinical routine with a proven high-end iterative algorithm

- Improve your process efficiency with an algorithm that can handle different metal implants
- Leverage your ability to address more challenging cases (e.g., dental fillings, pacemakers)
- Extend your referral base by delivering outstanding image quality in orthopedics

¹ iMAR is designed to yield images with a reduced level of metal artifacts compared to conventional reconstruction if the underlying CT data is distorted by metal being present in the scanned object. The exact amount of metal artifact reduction and the corresponding improvement in image quality achievable depends on a number of factors, including composition and size of the metal part within the object, the patient size, anatomical location and clinical practice. It is recommended, to perform iMAR reconstruction in addition to conventional reconstruction.

iMAR improves biliary drainage catheter visualization in ^{99m}Tc HIDA SPECT/CT following liver tumor embolization



- Patient with cholangiocarcinoma treated with right-sided dual-venous embolization
- Biliary drainage catheter draining externally left in situ
- Hepatobiliary scan with dynamic and static SPECT/CT performed to assess function of remnant liver tissue
- CT with iMAR eliminates metal artifact from biliary drainage catheter enabling clear visualization of catheter tip

Symbia Pro.specta™ SPECT/CT

SPECT

Scan acquisition:

60 stops per detector,
10 seconds/stop
OSEM3D 24i,4s 128x128 matrix

Injected dose:

^{99m}Tc -HIDA 3.8 mCi (142 MBq)

CT

Scan parameters:


120 kV/40ref mAs 0.33 sec rotation
32x0.7 mm collimation
2 mm recon
BR40 S3 iMAR

Data courtesy of Queen Elizabeth Hospital, Birmingham, United Kingdom.

iMAR is designed to yield images with a reduced level of metal artifacts compared to conventional reconstruction if the underlying CT data is distorted by metal being present in the scanned object. The exact amount of metal artifact reduction and the corresponding improvement in image quality achievable depends on a number of factors, including composition and size of the metal part within the object, the patient size, anatomical location and clinical practice. It is recommended, to perform iMAR reconstruction in addition to conventional reconstruction.


Symbia Pro.specta™ and its features are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Enhance your nuclear medicine department with intelligent oncology SPECT/CT imaging



CR Pro

Perform acquisitions with a continuous rotation to enable faster scan times.




AutoPlanar

Create virtual planar images from a single SPECT tomo acquisition; potentially reduce total imaging time.

Symbia Pro.specta SPECT/CT with myExam Companion CR Pro and Gated CR Pro: Faster scans with continuous rotation

Symbia Pro.specta™ is the first SPECT/CT to perform continuous-rotation acquisitions for all SPECT applications, including gated studies

- **CR Pro** reduces whole-body SPECT studies by 16 minutes¹
- **Gated CR Pro** reduces cardiac study time by up to 25%²



¹ Assumes a 5-bed study.

² Based on typical clinical workflow compared to a step-and-shoot gating acquisition.

Symbia Pro.specta and its features are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

Symbia Pro.specta SPECT/CT with myExam Companion AutoPlanar: automatically create planar images



SPECT/CT VRT



SPECT MIP



Virtual planar
anterior



Virtual planar
posterior

Automatically create supporting planar images without the need to perform a separate acquisition

- Up to 8 virtual static planar images (anterior, posterior, lateral, etc.) can be created from a SPECT acquisition
- Up to 2 virtual whole-body planar images (anterior and posterior) can be created from a whole-body SPECT acquisition

Děkuji za pozornost

Radomír Klíčník

Siemens Healthcare, s.r.o.

Sales & Portfolio Solutions Management - DI&AT&US

SHS EMEA CEECA CZE S&PSM-DI&AT&US

Podnikatelska 2924/2

612 00 Brno, CZECH REPUBLIC

Mobile: +420 602771300

<mailto:radomir.klicnik@siemens-healthineers.com>